

How to use DNA_Sequence_Detector

2020 August 31st

Shintaro Miyazaki
Satoshi Okada

- 1) The four files below have to be placed in the same folder.

SD_internal.ipynb

SD_external.ipynb

Sequence file 1 (query sequence, in FASTA format)

Sequence file 2 (subject sequence(s), in FASTA format, multi-FASTA format is acceptable)

- 2) The libraries below have to be installed.

numpy

tqdm

matplotlib

scipy

numba

datetime

- 3) Open the “SD_external.ipynb” file in Jupyter Notebook or Jupyter Lab.

- 4) Input the sequence file names. In the example below, the name of the query sequence file is “repeatUnit.txt” and the name of the subject sequence file is “testSequence.txt”.

```
repeatUnitFileName='CUP1_repeatUnit.txt'
f = open(repeatUnitFileName, 'r')
repeatUnitList = f.readlines()
f.close()
```

```
sequenceFileName='testSequence.txt'
f = open(sequenceFileName, 'r')
sequenceList = f.readlines()
f.close()
```

- 5) Run all cells.

- 6) DNA_Sequence_Detector returns 2 sets of output files.

SD.txt contains the result of detection of the query sequence. (in TSV format)

A number of PNG image files contain graphs showing the position of the query sequence in the subject sequence.